


Parallel evaporation of samples in 8 or 24 vials, tubes and microtitre plates

Features

- Precise control of inert gas flow combined with digital control of heating carefully evaporates your samples.
 - Interchangeable plates with either 8 or 24 hollow blowdown pins deliver an equal flow of gas to each tube, vial or well.
 - The absence of a vacuum avoids bumping, protecting the sample during evaporation.
 - Nickel-plated aluminium offers excellent chemical resistance.
 - Easy viewing of samples during evaporation.
 - Optional flowmeter precisely controls flow of inert drying gas.
- 
- Precise heat control
the absence of a
protects your s
and avoids bur

Control
and log your
chemistry



Enclosed design contains evaporating solvent, allowing subsequent trapping and collection of solvent via a high performance glass condenser

Can be converted into
GreenHouse Plus Parallel
Synthesiser

Typical Evaporation Times

Solvent	Samples	Volume	Vessel	Evap. Time
Methanol	24	2 ml	3.5 ml vial	20 mins
Methanol	8	5 ml	20 ml vial	22 mins
Acetonitrile	24	2 ml	3.5 ml vial	20 mins
Acetonitrile	24	2 ml	7 ml tube	35 mins
Water	8	5 ml	20 ml vial	157 mins
DMF	24	2 ml	3.5 ml vial	138 mins
DMF	24	2 ml	7 ml vial	145 mins

40°C Base Temperature. Flowrate 10 l/min (8 well), 20 l/min (24 well)

Compatible with:

- 7 ml GreenHouse tubes
- 8 or 24 position vial racks
- 13 mm, 13.8 mm, 15 mm, 24.3 mm & 27.8 mm \varnothing vials
- 8 or 24 well microtitre plates



Evaporate
8 vials, each
containing 5 ml of
methanol, in only
22 minutes

Safety relief valve on inlet prevents over-pressure during operation.

- Removable head features quick-release handles for easy exchange of Blowdown pin plates.

Select the appropriate Blowdown pin plate. Interchangeable plates with either 8 or 24 hollow pins deliver an equal flow of gas to each tube, vial or well.



Digitally controlled heating from the hotplate gently adds energy to the sample to speed the evaporation process.

Blowdown system with 24 pin plate, standard GreenHouse base, reaction block and 7 ml tubes

Insert the adapter into the GreenHouse base to accept vial racks or titre plates...

...or use the dedicated, low profile, Blowdown base.



Insert the adapter into the GreenHouse base



GreenHouse base and
24 position vial rack



GreenHouse base and
8 position vial rack



GreenHouse base and
24 position MTP



Blowdown base and
24 position vial rack



Blowdown base and
8 position vial rack

Blowdown base and
24 position MTP